




# DFMT FMT

Premier Radiographic Systems


www.AmRadMedical.com



Fast and Accurate Patient Positioning

### Digital Floor Mounted Tubestand

- Easy to read LCD shows SID continuous to table-top and table receptor with discreet SID of 40" and 72" to the wallstand
- Cross-table views are easily accomplished using the rotational and transverse positioning capabilities
- 10 ft. floor track allows for offtable and weight bearing exams and is customizable down to 8 ft.



Compatible with most DR Solutions



Convenient Lock-Release Pedals

# Standing Knees



Ergonomic Design

## **Elevating Four-way** Float Top Table

- Flat, four-way, floating top permits easy transfer and positioning of patients
- Designed for all types of applications, available in 400lb and 650lb patient capacities
- Four-pedal design allows for easy accessibility to table lock release
- Collision-protection electronics and fail-safe locks ensure safe operation

#### Options include:

- 35.5" wide table top
- Four-way lock release handle
- Non-elevating version available in 400lb and 650lb patient capacities

#### J1000 Wallstand

- Center mounting or side mounting of Bucky makes this wallstand flexible for any patient application or installation limitation
- Ergonomically designed hand control for ease of positioning wall receptor
- Side-mounted handgrips stabilize patients during PA exams

#### Options include:

- Overhead Grip
- Grid Cabinets in place of Bucky



35.5" Wide Table Top



Table Lock Release Handle



**Economical Alternative FMT** 



Overhead Grip for J1000 Wallstand

# System Options (L to R, above)

The AmRad™ Medical FMT and DFMT Systems can be configured to meet the needs of any department. Options include, but are not limited to:

- An extra-wide table top for the S222 elevating table to help patients feel more secure
- Four-way lock release handle for easy patient positioning can be located anywhere on the rails of the S223 non-elevating or S222 elevating tables
- The FMT is an economical alternative to the DFMT with the same versatility and functionality as the DFMT
- Patient Overhead Grip for use with the J1000 wallstand helps steady patients for upright exams



Call 1-800-729-9729 for additional information.



2901 W. Lawrence Avenue Chicago, IL 60625 USA

Phone: 773-353-4000 Fax: 773-588-3424 Toll Free: 1-800-729-9729

